FRED easy TRAINER

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Contents

1	General points	23
2	Intended use	25
3	Characteristics	25
4	Display and controls	26
4.1	Device	26
4.2	Electrodes and cable	28
5	Scenario selection procedure	29
5.1	List of scenarios	29
5.2	Selection of a scenario	29
5.2.1	Selection of a scenario on the device	29
5.2.2	Selection of a scenario with the infrared remote control	30
5.3	Description of the functions associated with keys on the infrared remote control	31
5.4	To switch the device off	31
6	Safety instructions	32
7	Cleaning and maintenance	32
8	Servicing and disposal at the end of the service life	32
9	Specifications	33
9.1	Physical and environment requirements	34
10	Warranty	34
11	Accessories	35
12	General algorithm	36
13	Spoken and displayed texts	37

WARNING!

This device does not deliver energies. Do not mistake it.

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1 General points

- The device complies with standard EN 55011 class B.
- The safety of the user and trouble free device operation shall only be guaranteed when original SCHILLER parts are correctly used. Further, only the accessories mentioned in these instructions for use or accessories that have been expressly authorised by SCHILLER may be used. If the device becomes damp as a result of temperature variations, it may only be put back into service once it is completely dry.
- The manufacturer shall not be held liable for the safety, reliability and performance of the device if:
 - Assembly, extensions, adjustments, modifications or repairs are not performed by the manufacturer or by personnel authorised by the manufacturer.
 - The electrical installation of the premises of use does not comply with locally applicable regulations.
 - The device is not used in accordance with the instructions for use.
- The manufacturer may supply, on request, circuit diagrams, lists of components, descriptions, calibration instructions or any other information, for use by the qualified technical personnel of the user for repairing such parts of the device as have been stated to be repairable by the manufacturer. Such supply shall not in any case constitute permission or approval to modify or repair the device.
- This document describes a device type that complies with the safety standards applicable at the time of printing. SCHILLER reserves all rights to the circuits, processes, names, software and devices mentioned in these instructions for use.
- The user guide is an integral part of the device and should always be kept near the device. Close observance of the information given in the user guide is a prerequisite for using the device as intended for and correct operation and ensures operator safety. Therefore, be sure to read the complete user guide.
- The user guide informs the device operator about the intended use, exact function, operation and required preventive maintenance. It is not a substitute for a product training.
- No part of these instructions for use may be reproduced without the written permission of SCHILLER.
- This device is not designed for any use other than that specifically described in this manual. Such use may be hazardous

Manufacturer

SCHILLER Médical SAS
4, rue Louis Pasteur – BP 90050
F-67162 WISSEMBOURG Cedex
Tel: +33 (0)3.88.63.36.00
Fax: +33 (0)3.88.94.12.82
E-mail: info@schiller.fr
Internet site: www.schiller-medical.com

Revision history:		
Version	Date	Comment
01-01	24.04.2003	1st issue
Août 2004	August 2004	1st revision
September 2004	September 2004	Addition of German version
March 2005	March 2005	Addition of alkaline battery
April 2005	April 2005	Modification of alkaline battery
May 2005	May 2005	Correction on side 20
May 2006	May 2006	New guidelines ERC / AHA 2005
October 2007	October 2007	Addition of AD folio
June 2008	June 2008	New remote control
April 2012	April 2012	Correction CE 0459 to CE and new symbols

2 Intended use

The FRED easy TRAINER is a simple cost-effective simulator of the FRED easy that is designed to be used exclusively for training. In no event may it be used as a conventional defibrillator. It should not be connected to a patient. It should only be used on a mannequin. It has been designed for the purpose of training users of SCHILLER's FRED easy automated defibrillator. The device provides several simulations or scenarios (from 0 to 8) that aim to help to familiarise with the FRED easy and shows that the users have the knowledge required to use it in case of emergency. The device also allows the trainer to work using an infrared remote control to select a scenario or to intervene in real time on the course of an intervention.

3 Characteristics

The FRED easy TRAINER allows the **simulated** delivery of electric shocks. There is no high voltage in the device, and training is risk free.

The FRED easy TRAINER may be used on all types of mannequins as it does not deliver an electric shock. It may only be used with the training electrodes sold by Schiller.

The FRED easy TRAINER may be used with an infrared remote control. The remote control allows the trainer to control the training actions.


The FRED easy TRAINER has 9 pre-recorded training scenarios which simulate realistic cardiac arrest configurations.

The basic equipment of the FRED easy TRAINER consists of :

- 1 rechargeable 6 Volts battery and a battery charger or 1 alkaline 6 Volts battery
- Training electrodes consisting of a single use set of 2 electrodes and a reusable cable. The connector of the cable is fitted with a magnet allowing the device to recognise the presence of the electrodes in the device without any electrical contact
- 1 remote control
Note: The remote control is programmed at delivery and ready for use.
- 1 AD (automatic defibrillation) repositionable folio.

4 Display and controls

4.1 Device

- 1 Green On / Off button and for selection and save of the configuration
- 2 Battery housing
- 3 Remote control ; only the keys 0 to 9 are used
- 4 LCD display : displaying of elapsed time, number of shocks simulated, battery autonomy, text messages
- 5 Orange button for shock simulation according to the configuration of the device
- 6 Infrared receiver
- 7 Green diode for correct operation simulation
- 8 Orange diode for electrode fault simulation and scenario selection by remote control
- 9 Connector for connecting cable of the training electrodes
- 15 Infrared transmitter

Symbols used on the device


Symbol for recognition of electrical and electronic equipment.
The device must be disposed of in a municipality approved collection point or recycling centre when it is no longer required.
Improper disposal harms the environment and human health due to the presence of dangerous substances in electrical and electronic equipment.
Manufacturer symbol, manufacturing date

Symbols used on the battery

	The unit/component can be recycled.
	Battery may not be disposed of with domestic refuse.
	Do not burn, saw up or crash the battery
	Rechargeable battery
	Do not short the battery
	Storage temperature for the battery: Unlimited; 0...+50 °C
	Expiration date

4.2 Electrodes and cable

- 10 Adhesive training electrodes; single use
CAUTION! : The electrodes are not designed to connect a patient electrically to the device; they are isolated. Don't use them with the defibrillator FRED easy!
- 11 Connector of the cable with a magnet to recognise the electrodes
- 12 Connection cable between electrodes and device
CAUTION! : The cable doesn't be used with the defibrillator FRED easy!
- 13 Discs of the cable to stick on the electrodes
- 14 Self adhesive discs of the electrodes

5 Scenario selection procedure

5.1 List of scenarios

SNR Shock not recommended
SR Shock recommended
Elect prob Electrode problem

Scenario 0	SNR
Scenario 1	SNR / SR
Scenario 2	SNR / SR / SNR
Scenario 3	SR
Scenario 4	SR / SNR
Scenario 5	SR / SNR / SR
Scenario 6	Elect prob / SR
Scenario 7	SR/ Elect prob / SR
Scenario 8	SR / SR / SNR
Scenario 9	Infrared control function

Note : the last step of each scenario is indefinitely repeated.

5.2 Selection of a scenario

5.2.1 Selection of a scenario on the device

- Press the green button (1) to switch the device on. It functions with the last selected scenario.
- Place the device on standby mode by pressing the green button (1) for a few seconds. The green diode (7) flashes.
- Press at first the orange button (5), keep it pressed and press then the green button (1) to access in the configuration menu.

The configuration menu is displayed on 2 successive screens of 4 lines.

The choice in a row is done using the green button (1) and the transition to the next line is done with the orange key (5).

- Line 1.1: choice between the semi-automatic operation mode and the automatic operation mode of the device.
Note: in automatic operation mode, position on the front plate the folio delivered with the device.
- Line 1.2: this line appears only in semi-automatic mode: choice between automatic or manual analysis start.
- Line 1.3: You can go to the next page of the menu.
- Line 1.4: Saving the configuration and exit the menu (with the green button).

- Line 2.1: Activating the metronome function
- Line 2.2: this line appears only when the metronome function is activated: adjusting the frequency of the metronome.
The frequency of the metronome can be set between 80 and 150 p/min by steps of 5 p/min.
- Line 2.3: choice among the 9 scenarios (0 to 8) or infrared function (scenario 9).
- Line 2.4: Saving the configuration and exit the menu (with the green button).
- Select in this way the scenario according to the configuration of the FRED easy which you wish to use for training. Validate your selection by pressing the green button (1). The device then switches off and goes in standby mode. The green diode (7) flashes to inform the user that the device is operating normally.

The device can then be switched on normally and it will operate according to the selected scenario.

5.2.2 Selection of a scenario with the infrared remote control

The trainer has the possibility to select a scenario with the remote control (by pressing the button on the remote control corresponding to the selected scenario) when the device is in the standby mode (device off and green diode is flashing).

The orange diode (8) under the electrode connector flashes in the moment of the scenario selection to inform the user that the scenario has been taken in account by the device.

The device can then be switched on normally and it will operate according to the selected scenario.

In the case of the selection of scenario number 9, the device waits for the remote control's orders to change its status. If the device doesn't receive orders from the remote control, the shock will always be recommended like in scenario 3.

5.3 Description of the functions associated with keys on the infrared remote control

Key 0	Immediate move to the next step
Key 1	Restart at the analysis step; reset of the counters of time and number of delivered shocks
Key 2	Shock recommended
Key 3	Shock not recommended
Key 4	The rhythm has changed
Key 5	Press on the electrodes (bad contact)
Key 6	Movement detected. Stand clear.
Key 7	Simulation of device off and green diode (7) is flashing CAUTION! The device is now in the standby mode. Press on the green button (1) to switch the device on.
Key 8	Resume CPR
Key 9	Pause: the current scenario is interrupted: the display is frozen and the green diode (7) does not flash. Press again on key 9 to continue the training

Note: The keys 2 to 6 should be activated before the end of the analysis phase to be taken in account.

5.4 To switch the device off

Press at the same time for a few seconds the orange button (5) and the green button (1), then slacken them to switch the device completely off. It is then no longer consuming energy.

CAUTION:

The device is only completely switched off when the green diode is no longer flashing. Do not leave it in standby mode, as this empties the battery unnecessarily.

6 Safety instructions

The following points shall be specially respected:

- The device has been designed for training purposes **only**. It may in any case not be used for defibrillation.
- The device is used with training electrodes.
- These instructions for use describe the complete FRED easy TRAINER device with all its optional features.
- The device is not designed for any use other than that described in these instructions. Such use is forbidden.
- The user shall be entirely responsible for use of other accessories than from SCHILLER.
- The user must be familiar with the functioning of the device.
- Before each use, the user must ensure the device is working correctly.
- The device is not designed for use in areas where an explosion hazard may occurs

7 Cleaning and maintenance

The outer surfaces of the device may be cleaned with cotton wool moistened with alcohol.

CAUTION!

**Risk of damage to the device.
Never use phenol based cleaning agents or agents containing
peroxide derivatives to disinfect the surfaces of the device housing**

During cleaning, the batteries must be removed from the device. No liquid may be allowed to enter the device. If liquid does penetrate into the device, clean the device and have it verified completely.

8 Servicing and disposal at the end of the service life

Check the device performance once a year.

At the end of its service life, the device and the accessories must be disposed of in compliance with the local regulations.

It is mandatory to recycle the parts of the device separately and return them to the available collection and recycling centres.

Apart from the plug-in battery, the device does not contain hazardous material and can be disposed of like any other electronic equipment.

9 Specifications

- Operating mode: training in semi-automatic and automatic defibrillation
- Power supply:

Alkaline battery 6 Volts composed of 4 alkaline batteries LR6 1,5 V

CAUTION!
Do not recharge! Do not put in the charger!
Do not use the battery in the defibrillator FRED easy!

Rechargeable Ni /Cd battery 6 Volts, 650 mAh

CAUTION!
Do not use the battery in the defibrillator FRED easy!

ATTENTION!
RISK OF DETERIORATION OF THE DEVICE:
Do not plug the 12V battery of the defibrillator FRED easy in the training device FRED easy TRAINER!

Battery charger 100-240 VAC,
50-60 Hz, 0,35A max ; Class II ;
Charge time of an empty battery :
1h

Battery adapter

1 to 100% on-screen battery indicator corresponding to the battery voltage range from 4.5 Volts minimum to 6 Volts.

- Power consumption of the device: 35 mA without voice; approximately 250 mA with voice.
- Detection of electrode presence by means of a magnet placed in the connector of the connecting cable for the training electrodes.

9.1 Physical and environment requirements

- Transport /storage:
Temperature: – 30...+ 50°C
Relative humidity rate: 30...95 %, without condensation
Atmospheric pressure: 500...1060 hPa
- Operation:
Temperature: 0...+ 50 °C
Relative humidity rate: 30...95 %, without condensation
Atmospheric pressure: 700...1060 hPa

Dimensions and weight

Width	230 mm
Depth	220 mm
Height	70 mm
Weight	< 1 kg

10 Warranty

SCHILLER guarantees its devices to the purchaser against all material or manufacturing defects for a period of one year as from the original purchase date. During this one-year period, SCHILLER, as soon as it receives a device presenting a material or manufacturing defect returned accompanied by written notification of the defect, will repair or replace the defective parts or the complete product if necessary.

All postal, shipping or handling costs will be borne exclusively by the user.

The manufacturer shall supply, on request, circuit diagrams, lists of components, descriptions, calibration instructions or any other information, for use by the qualified technical personnel of the user for repairing such parts of the device as have been stated to be repairable by the manufacturer. Such supply shall not in any case constitute permission or approval to modify or repair the devices.


All rights are reserved for the circuits, processes and names mentioned in these instructions for use.

SCHILLER shall not be held liable as a result of this warranty for any secondary or indirect damage, or in the event of unauthorised repairs or modifications having been made or undertaken, or when the whole or a part of the product has been damaged by accident, incorrect use or deliberate misuse. This warranty does not cover the battery, normal wear and tear, marks or any other irregularity of appearance which does not prevent or hinder the performance of the device.

11 Accessories

Part Number	Designation
6-90-0053	Case for 4 alkaline batteries (without batteries)
53108	Set of 4 alkaline batteries 1,5 Volt LR6
4-07-0010	Rechargeable NiCd 6 Volts battery
3-55-0030	Charger for 6 Volts NiCd battery, 100-240 VAC, 50-60 Hz Europe
0-21-0006	Reusable cable for training electrodes 0-21-0005. Length: 2 m.
0-21-0005	Adhesive training electrodes (10 sets of 2 electrodes), single use.
0-99-0004	Remote control, ready for use.
9-73-0000	Set of 4 alkaline batteries 1,5 Volt LR3 for the remote control.
6-06-0012	Repositionable Automatic Defibrillation folio.

12 General algorithm

Note :

⏮ 1 ↻ 20 seconds = Text 1 repeated every 20 seconds

Rem 1: The step "PRESS ANALYSE" (text 2) is cancelled if the analyse is automatic or if the device functions in the automatic mode.

Rem 2: The texts 14, 15 and 16 are only present when the device operates in the automatic mode.

13 Spoken and displayed texts

	Spoken text	Displayed text
1	Place electrodes on chest and plug into machine.	PLACE ON ELECTRODES PLUG INTO MACHINE
2	Press green button.	PRESS GREEN BUTTON
3	Do not touch the patient, analysing.	DO NOT TOUCH THE PATIENT ANALYSING
4	Do not touch the patient, charging.	DO NOT TOUCH THE PATIENT CHARGING
5	Stand clear of patient, press orange button.	PRESS ORANGE BUTTON TO SHOCK
6	No shock advised	NO SHOCK ADVISED
7	Heart rhythm has changed.	RHYTHM HAS CHANGED
8	Immediately resume CPR. 30 chest compressions then 2 rescue breaths; continue until patient is breathing normally.	30 CHEST COMPRESSIONS THEN 2 RESCUE BREATHS
9	Movement detected. Stand clear.	MOVEMENT DETECTED STAND CLEAR
10	Check the electrodes.	CHECK THE ELECTRODES
14	Warning! Shock will be delivered 5-4-3-2-1.	WARNING! SHOCK
15	Make sure the patient is unresponsive. If not, switch unit off.	MAKE SURE THE PATIENT IS UNRESPONSIVE. IF NOT SWITCH UNIT OFF.
16	Do not touch the patient, analysis will begin	DO NOT TOUCH THE PATIENT ANALYSIS WILL BEGIN

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